

# nüvi cameras

every photon counts

## h·nüTDI

HIGH SENSITIVITY  
TIME-DELAY  
INTEGRATION  
MODE

### A NEW STANDARD FOR LOW LIGHT IMAGING

### NÜVÜ™ TDI CAMERA OPTIMIZED FOR DEMANDING SCANNING APPLICATIONS



#### OUTSTANDING LOW-LIGHT SCANNING PERFORMANCES:

Line rate up to 103 kHz  
Resolution of 4096 (H) x 128 (V) pixels  
Average readout noise lower than 65 electrons


#### TIME-DELAY INTEGRATION (TDI) PRINCIPLE:

A readout mode based on the concept of the accumulation of cumulative exposures of the same object as it is moving linearly under the detector. It synchronizes the transfer of the charges from one line to the next with the same speed as the object is moving under the camera.



**Fig 1.**  
Example of imaging the word NÜVÜ using TDI readout mode. The intensity of the signal is increased as the word moves across the detector

#### TYPICAL SPECTRAL RESPONSE



**Fig 2.**  
Typical spectral response as a function of wavelength, as specified by the detector manufacturer


# SIMPLE INTEGRATION & QUALITY PRIORITY

Nüvü Camēras offers the highest standard of technology in a compact thermoelectrically cooled camera. The technology at the heart of the HNü was originally designed for space exploration where the need for state-of-the-art instruments drives innovation. Now optimized and extended to a broad range of applications, the user- friendly HNü provides many advantages to efficiently bridge the gaps between purchase, setup, discoveries and usage.

All parts are treated in compliance with high vacuum requirements, including all metal sealed in a Class 10,000 cleanroom to ensure the longest camera lifetime without maintenance. Nüvü Camēras uses at least  $\lambda/10$  quality windows, essential for optimal image quality.

Consultation services are available on demand.

## TECHNICAL DRAWINGS



# h·nü TDI

| CHARACTERISTICS                 | SPECIFICATIONS                 |
|---------------------------------|--------------------------------|
| Line Rate                       | 103 kHz                        |
| TDI Transfer Direction          | Unidirectional                 |
| Interface                       | Camera Link Extended Full      |
| Resolution                      | 4096 x 128 pixels              |
| Data Rate                       | 30 MHz per output              |
| Pixel Size                      | 12 x 12 µm                     |
| Data Format                     | 14 bits                        |
| Outputs                         | 16                             |
| Size (H W D)                    | 155.6 x 155.6 x 177.8 mm       |
| Mass                            | 5 kg                           |
| Operating Temperatures          | 0°C to 30°C                    |
| Power Supply                    | 100-240 V, 50-60 Hz, max. 3 A  |
| QE                              | See spectral graph (Fig.2)     |
| Sensor Type                     | Back thinned                   |
| Effective Area                  | 49.152 mm x 1.52 mm            |
| Readout Mode                    | Conventional                   |
| TDI Line Rate Control           | Internal or external           |
| Full Well Capacity <sup>1</sup> | 70k electrons                  |
| Readout Noise <sup>1</sup>      | 65 electrons                   |
| Binning                         | Vertical, 1 to full binning    |
| Image Processing                | None                           |
| Sensor Cooling                  | Thermoelectric, to sub-zero °C |
| Dark Current <sup>1</sup>       | <130 ē/pixel/s @ 0°C           |

1 Typical values, these numbers may vary depending on the CCD detector.



Contact us at:  
info@nucameras.com  
+1 514 733 8666  
Montreal (Quebec)  
CANADA

HNü and NüPixel are the intellectual property of Nüvü Camēras. All other brands are properties of their respective owners. Incremental changes are made to the products and specifications are subject to modification without prior notice.

HNü TDI Specification Sheet 1.2.1  
© Nüvü Camēras, 2022